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BALANSSEN PRODUCT CARE AND CLEANING

Please see the following information so that you can enjoy the premium properties of Balanssen medical products as long as possible:

- Balanssen products are filled with organic husk of Buckwheat, Spelt or Millet that should not have contact with water.
- Pillow/roller covers or cases are made of natural fabric. You can wash them according to the instruction on the label and solely after removing the filling. We recommend pouring the husk into a container before washing.
- You do not clean the seed husk. However, we recommend airing it out recurrently which extends the healthy values of the filling.
- In case a Balanssen product contacts water one should remove the seed husk and dry it at an airy and warm place.

USING BALANSSEN PRODUCTS

- All Balanssen medical products have been filled with seed husk in an optimal way. Considering individual preferences of patients as for the rigidity and the height of pillows, Balanssen products are equipped with zippers that enable the adjustment of the filling amount.
- The adjustment of the filling amount makes it also possible to use the same product in case of a new indication, e.g. a roller might be used for sleeping, rehabilitative exercises, elevation of the lower limbs or unweighting the lumbar spine.
- The seed husk that shifts around adapts to the shape of a user's body and supports it stably. In case of the sleeping pillows, before lying down, we recommend initial forming of a dimple under the head and forming a roller under the neck which enables using the same product whether one sleeps on their side, back or stomach.