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BALANSSEN PRODUCT CARE AND CLEANING

Please see the following information so that you can enjoy the premium properties of Balanssen medical products as long as possible:

- Balanssen products are filled with organic husk of Buckwheat, Spelt or Millet that should not have contact with water.
- Pillow/roller covers or cases are made of natural fabric. You can wash them according to the instruction on the label and solely after removing the filling. We recommend pouring the husk into a container before washing.
- You do not clean the seed husk. However, we recommend airing it out recurrently which extends the healthy values of the filling.
- In case a Balanssen product contacts water one should remove the seed husk and dry it at an airy and warm place.

USING BALANSSEN PRODUCTS

- All Balanssen medical products have been filled with seed husk in an optimal way. Considering individual preferences of patients as for the rigidity and the height of pillows, Balanssen products are equipped with zippers that enable the adjustment of the filling amount.
- Excessive point pressure exerted on body soft tissues inhibits blood flow and leads to tissue death, which, as a consequence, may result in bedsores and chafe.

Purpose of the pillows and the seats is relieving the places that are exposed to long-lasting pressure through natural tipping of the husk filling and letting the air flow easily between patient's body and a product, thus preventing the incidence of bedsores and chafe.

The products are recommended for patients who spend much time seating or they are required to lie for a sustained period of time.